

# Mass Isotopologue Distribution Analysis (MIDA)



Van Andel Research Institute Metabolomics Core  
Brejnev Muhire, Christine Isaguirre and Ryan Sheldon

MIDA is a computer program used for isotope natural abundance correction using either labeled samples or theoretical mass isotopologues distribution values (MIDs). It generates natural abundance corrected values in comma-separated format and high-quality figures. It is written in R and is equipped with a user-friendly ShinyApp graphical interface. ([pub.link](#))

## 1. Program structure

- app.R, server.R, ui.R : source code graphical user-interface
- mida.R : R script used to run MIDA in command line mode
- utils.R : functions used throughout the program
- MIDA.bat : batch file used to launch the MIDA ShinyApp on Windows OS

## 2. Input format

Input data consist of table of compound concentration values and a metadata describing samples and file listing names of all compounds in .csv, xls or xlsx format.

Name	Data File	Type	Level	Acq. Date-Time	Pyruvate Results				Lactate Results			
					Resp.	Pyruvate M+1 Results	Pyruvate M+2 Results	Pyruvate M+3 Results	Resp.	Lactate M+1 Results	Lactate M+2 Results	Lactate M+3 Results
Blank_1	Blank_1.D	Sample		7/15/2019 16:31	15380	1980	614	1216	163959	32336	377893	40380
0min_a	0min_a.D	Sample		7/15/2019 17:14	6442646	868006	324565	24364	5217190	1250544	830266	103262
0min_b	0min_b.D	Sample		7/15/2019 17:57	5020177	662908	240234	16755	5064552	1234812	821446	119660
0min_c	0min_c.D	Sample		7/15/2019 18:40	3717213	478084	174101	12140	3601465	849116	698383	75078
Blank_2	Blank_2.D	Sample		7/15/2019 19:23	29424	2561	928	1353	333694	66642	416185	18587
20min_a	20min_a.D	Sample		7/15/2019 20:06	6701832	936370	344214	26422	4922736	1189268	821469	128954
20min_b	20min_b.D	Sample		7/15/2019 20:49	8707993	1223262	457281	37342	5477819	1295587	858876	99288
20min_c	20min_c.D	Sample		7/15/2019 21:32	7724890	1075413	388054	32863	4630630	1075568	797195	84972
Blank_3	Blank_3.D	Sample		7/15/2019 22:15	35560	4736	1959	76	394830	76273	420337	28706
1h_a	1h_a.D	Sample		7/15/2019 22:58	5181952	685477	258728	22756	4536691	1088041	797703	100957
1h_b	1h_b.D	Sample		7/15/2019 23:41	7096782	953857	363127	31163	5197828	1226366	842972	
1h_c	1h_c.D	Sample		7/16/2019 0:24	7841484	1044085	398633	33100	5483322	1347829	879188	201018
Blank_4	Blank_4.D	Sample		7/16/2019 1:07	32889	3144	1115	1482	487766	109553	414711	15110
2h_a	2h_a.D	Sample		7/16/2019 1:50	9105700	1303694	483999	43042	5864504	1409711	906430	116367
2h_b	2h_b.D	Sample		7/16/2019 2:33	6556746	873036	323530	30224	5589634	1333775	877744	105710
2h_c	2h_c.D	Sample		7/16/2019 3:16	7121191	951436	365710	30754	5431881	1359129	857430	137872

Fig. 1: Uncorrected data

Name	Labeling
Blank_1	blank1
0min_a	unlabeled
0min_b	unlabeled
0min_c	unlabeled
Blank_2	blank2
20min_a	labeled_sample1
20min_b	labeled_sample1
20min_c	labeled_sample1
Blank_3	blank3
1h_a	labeled_sample2
1h_b	labeled_sample2
1h_c	labeled_sample2
Blank_4	blank4
2h_a	labeled_sample3
2h_b	labeled_sample3
2h_c	labeled_sample3

Fig. 1: Meta Metadata

Compound	Included
Pyruvate	1
Lactate	1
L.Alanine	1
L.Glycine	1
X3.Hydroxybutyrate	0
L.2.Aminobutanoic.Acid	0
Beta.Alanine	0
N.Acetyl.glutamate	0
L.Valine	0
L.Leucine	0
L.Isoleucine	0
Succinate	1
L.Proline	1
Fumarate	0
Pyroglutamate	0
L.Methionine	1
Serine	1
alpha.Ketoglutarate	0
L.Threonine	0
L.Phenylalanine	0
Malate	1

Fig. 3: Compound list

### Note:

All files need to in the same input directory and have the exact matching prefix but different suffix (“\_uncorr.csv”, mtd.csv, and ‘cmpd.csv’)

Eg: “BHB\_uncorr.csv”, “BHB\_mtd.csv”, “BHB\_cmpd.csv”

Multiple experiment file can be put in one directory given that the have different prefix name.

In compound list 1 and 0 respectively indicate compounds which should or should not be corrected.

...

### 3. Workflow

- **Running MIDA**

Place the input \*\_uncorr.csv files and the associated \*\_mtd.csv and \*\_compd.csv files in one folder. Within the MIDA App GUI use the “Choose folder” button to select the input folder.

The settings tab allows to set the following:

- Correction method: whether to use labeled samples or theoretical MIDs
- Baseline detection level: The threshold detection level at which compounds are filtered
- List of compounds: correct for all compounds or a user-selected list (in “\_compd.csv”)
- Running mode: skip or overwrite the existing output directories of processed input files

To run the command using the command prompt

```
Rscript mida.R InputPath List 2000 overwrite or  
Rscript mida.R InputPath All 1000 skip
```

- **Data processing**

1. Reads / checks input table for samples, metadata and compound list
  - Check the ordering of isotopes (M+0, M+1, M+2...M+n)
  - Select compounds as indicated in the compound list “\_compd.csv” file
  - Read metadata and check whether samples names match those in the input table.
2. Substitute missing values with zero
3. Imputation of missing values (per group, per isotopologue, min val > 0)
4. Remove all compounds below baseline detection (eg. 2/3 of unlabeled samples M+0 are below 1000 or all labeled sample entries are < 1000)
5. Natural abundance correction: for each compound, create the correction matrix CM, from unlabeled data (for multiple sample, the mean of each isotopologue is used). Each labeled sample is corrected using the following formula:

$$\mathbf{A}_{\text{corr}} = [\mathbf{CM}]^{-1} \times \mathbf{A}_{\text{obs}}$$

$[\mathbf{CM}]^{-1}$ : inverse of correction matrix,  $\mathbf{A}_{\text{obs}}$ : labeled sample MIDs,  $\mathbf{A}_{\text{corr}}$ : corrected MIDs

Note: CM can also be obtained from theoretical MIDs for a given compound

6. Removing compounds for which 2/3 or more of unlabeled samples (M+0) have total fraction enrichment <95 or >105.
7. Generating figures and tables.

### 4. Output

The output is generated in same folder as input file with name prefix matching corresponding input plus a suffix “\_out”. Eg. “BHB\_out”. The output consists of the following:

**Figures:**

Total Enrichment Bar Plots.  
Percent Enrichment Bar Plot  
TCA Cycle show a bar plot of each compound  
Heatmap compound x samples  
Coefficient of variation across samples

**Tables:**

Corrected MIDs in .csv format  
Excel formatted data for input for PRISM

The figures and tables can be visualized using the program interface. Example of figures are shown below:




Fig 4 Total enrichment barplot




Fig 5 Percentage enrichment barplot.




Fig 6 TCA cycle plot.




Fig. 7: heatmap




Fig. 8: Input total pool size coefficient of variation